

# MediaRecorder 6.5 - Quick Start Guide

- [Welcome to MediaRecorder](#)
- [Install MediaRecorder](#)
- [Set up MediaRecorder](#)
- [Options](#)
- [Record video](#)
- [For more information](#)

# Welcome to MediaRecorder

## MediaRecorder 6.5

Information in this document is subject to change without notice and does not represent a commitment on behalf of Noldus Information Technology bv. The software described in this document is furnished under a license agreement. The software may be used or copied only in accordance with the terms of the agreement.

Copyright © 2023 Noldus Information Technology bv. All rights reserved. No part of this publication may be reproduced, transmitted, transcribed, stored in a retrieval system, or translated into any other language in whole or in part, in any form or by any means, without the written permission of Noldus Information Technology bv. MediaRecorder is a product of Noldus Information Technology bv. Other product names are trademarks or registered trademarks of their respective companies.

Documentation: Olga Krips, Viktor Kovacs, Leanne Loijens  
December 2023

Noldus Information Technology bv.  
International headquarters  
Wageningen, The Netherlands  
Telephone: +31-317-473300  
Fax: +31-317-424496

E-mail: [info@noldus.nl](mailto:info@noldus.nl)

For addresses of our other offices and support, please see our website  
[www.noldus.com](http://www.noldus.com)

# Install MediaRecorder

**IMPORTANT** - If you have an older version of MediaRecorder - Uninstall MediaRecorder 6.0 and older versions before you install MediaRecorder 6.5.

## To install MediaRecorder 6.5

1. Download the installation files from [my.noldus.com](http://my.noldus.com), and extract the contents.
2. Double-click the file **MediaRecorder [Version number] Setup.exe**.
3. Under Installation type, select **Standard**.
4. In the **Drivers and tools** field, choose the driver of your video device.
5. Click **Next**, accept the End-User License Agreement and click Install.

### Notes

- See [Specifications](#) for the supported devices.
- See [Recording Devices](#) for specific information on your device like the maximum number of cameras you can use simultaneously and the maximum supported frame rate and resolution.

# Set up MediaRecorder


## up and running quickly

This Quick Start Guide guides you through the main steps to record video with MediaRecorder. This Quick Start Guide contains basic information on MediaRecorder functionalities. It does not contain information how to connect and setup your cameras. Also, it does not contain information on the prerequisites and restrictions for specific camera types. You can find this additional information in the MediaRecorder Help. Press F1 in the program to open it, or type MediaRecorder Help in the Windows apps screen.

## start mediarecorder


1. Connect all cameras and microphones to the computer with MediaRecorder.
2. Depending on your license, do one of the following:
  - If you have a hardware key - Insert it into a USB port of the computer and start MediaRecorder.
  - Start MediaRecorder and activate your software license key. You can choose between a Floating or a Fixed activation.
    - **Floating** — Choose Floating if you want to be flexible on which computer you use MediaRecorder and your computer is connected to the internet.
    - **Fixed** — A Fixed activation is linked to one computer. Your computer must have an internet connection to activate/deactivate a fixed license. After you have activated the license you can use MediaRecorder without internet.

The first time MediaRecorder opens, the preview window is black. That is because you still have to select the cameras. Proceed with the next steps to do so.




## select cameras

1. Choose **File > Video Settings**, or click **Video Settings** in the main window.



**TIP** If the **Video Settings** button is not visible, the recording properties may be closed. Click **Recording properties** at the bottom of your window.





2. Select your camera from the **Video device** list.



	Use	Video name	Video device	F
1	<input checked="" type="checkbox"/>	Video 1	Select video device...	

See [Set up IP cameras in the MediaRecorder](#) Help how to select IP cameras.

3. Select the frame rate and resolution.

	Use	Video name	Video device		Frame rate	Resolution	Audio
1	<input checked="" type="checkbox"/>	Video 1	Microsoft LifeCam Studio		30.00	1280 x 720	Audio
2	<input type="checkbox"/>	Video 2	Select video device...				
3	<input type="checkbox"/>	Video 3	Select video device...				



4. For advanced camera settings, click the video symbol next to the name of your camera.  
The available options depend on your camera. For some cameras you can zoom and adjust settings for brightness and contrast.

	Use	Video name	Video device		Frame rate
1	<input checked="" type="checkbox"/>	Video 1	Microsoft LifeCam Studio		30.00
2	<input type="checkbox"/>	Video 2	Select video device...		

5. Optionally change the names in the **Video name** column, to, for example, Overview camera, or Close-up participant 1.
6. To record audio as well, continue with [select microphones](#).
7. Optionally, select some additional settings for your video files. See [Options](#).
8. When satisfied with the settings, click **OK** to leave the **Video Settings** window and start recording. See [Record video](#).

## select microphones

1. Select the microphones from the **Audio device** list. Choose between your computer microphone, a separate microphone, or the audio stream of your camera. Choose **No audio** to not record audio.

	Audio name	Audio device	
	Audio 1	No Audio	
	Audio 2	Headset Mic (Realtek High Definition	
	Audio 3	Jack Mic (Realtek High Definition Au	
	Audio 4	Microphone Array (Realtek High Def	
	Audio 5	RTSP IP Camera (Noldus RTSP Source	
		No Audio	

2. Optionally, record the same audio stream in all videos. To do so, select this checkbox. (Please note that in Picture-in-Picture mode, you may select which audio source to use for all the videos. See [Select an audio device](#)).

### Output settings

Output mode: Separate videos

PiP Position: Top - Left

PiP Size: 20%

PiP Margin: 5%

1

2

3


4

☒ Use audio device 1 for all videos

3. Optionally, click the microphone button for advanced audio settings.

Settings

Frame rate	Resolution	Audio name	Audio device
5.00	1280 x 800	Audio 1	Microphone (High Definitio...
0.00	1280 x 800		


 Advanced audio settings

The available options depend on your microphones. Move the sliders to adjust the recording volume.

4. Optionally, change the names in the **Audio name** column.

## To check the recording volume

1. Right-click the speaker icon in the bottom-right corner of your window and select **Recording devices**.



2. Check that the green bar is maximal when speech level is normal.





3. Continue with [Options](#), or click **OK** to close the Video Settings window and [start recording](#).

#### Note


Audio is needed and has been tested only for cameras that are supported with The Observer XT. See **Recording Devices** in the MediaRecorder Help which cameras are.

# Options

For a full overview of the options, see the MediaRecorder Help.

## set a recording time

1. Choose **File > Video Settings**.
2. Select **Use maximum recording time and enter the time**.



The screenshot shows the 'Timer settings' dialog box in MediaRecorder. It has a sidebar on the left with a dark grey bar. The main area contains the following settings:

- Time format:** A dropdown menu showing 'HH:mm:ss'.
- Decimals:** A dropdown menu showing '0'.
- Example:** A text field showing '16:43:52'.
- Use maximum recording time:** A checkbox that is checked, highlighted with an orange border, with a text field next to it showing '00:30:00'.

At the bottom right, there are three buttons: 'OK', 'Cancel', and 'Apply'.

## select output settings for multiple cameras

If you have selected more than one video source, MediaRecorder can create a separate video file for each camera, but also one video of all camera streams. Choose **File > Video Settings** and select one of the options under **Output Settings**. The following options are available:

- **Separate videos** — One video file for each camera.
- **Picture by Picture** — One video file with the images from each camera next or above each other.



- **Picture in Picture** — One video file with a large image of the first camera with the images from the other cameras embedded in it.



If you select this option the following settings become available:

- **PiP Position** — The position of the embedded images.
- **Pip Size** — The size of the embedded images. The percentage is based on the width of the first camera image. The aspect ratios of the embedded camera images are maintained.
- **Pip Margin** — The margin of the embedded images. When the percentage is zero, the embedded images are positioned directly on the edge of the main image. The percentage is based on the width of the first camera image.

#### Note

If possible, choose the same frame rate for all the videos for Picture in Picture or Picture by Picture. If you choose different frame rates the videos with the lowest frame rates may flicker.


# manage video files

## To select the video file location

The default location where MediaRecorder stores its video files is:

- D:\Noldus\MediaRecorder\Video Files — If Noldus IT installed your system
- C:\Users\Public\Public Documents\Noldus\MediaRecorder\Video Files — In all other cases.

To change the default location, click the button next to the **File location** field of the main MediaRecorder window and select another folder.




You must have writing rights in this folder. Also make sure that there is enough disc space. Dependent on the recording device, frame rate, resolution and amount of movement, MediaRecorder videos can be up to 600 Mb per hour.

## To select a file base name

The name of the recorded video file consists of the **File base name** followed by:

- The video input name.
- The date and the time the videos were recorded. The format of the date and time depends on the computer's regional settings.
- The number of the video camera.



## save the settings

The settings you create for the different cameras are automatically saved when you close MediaRecorder. They are saved to the file **Current settings.mrs** that is present in the folder C:\Users\Public\Public Documents\Noldus\MediaRecorder. These settings are used when you restart MediaRecorder. But you can also save your own settings, or save different settings for different setups, for example, a webcam for a test participant's face in a usability test room and three overview IP cameras in an observation room. This way you do not have to select all cameras manually each time.

## Procedure

1. Choose **File > Save Configuration As**.
2. Give the configuration a name and click **Save**.

## To open a saved configuration


Choose **File > Open configuration** and select a **.mrs** file.

### Notes

- **important** Never shut down the computer by pressing the power button or cutting off the electricity. Always use the Windows shut down feature to close your computer.
- To go back to the default configuration, choose **File > Open Configuration**. Open the file **Default Settings.mrs**.

# Record video

When you selected all cameras, you see a preview in the main MediaRecorder window.



**TIP** If you selected multiple microphones, it can be useful to mute all but one audio channel in the preview. Select the audio source in the **Mute** column of the table under the preview. Muting audio in the preview does not influence the recording.

Nr	Use	Video name	Audio name	Mute	Volume
1	<input checked="" type="checkbox"/>	Video 1	Audio 1	<input type="checkbox"/>	
2	<input checked="" type="checkbox"/>	Video 2	Audio 2	<input checked="" type="checkbox"/>	
3	<input checked="" type="checkbox"/>	Video 3	Audio 3	<input checked="" type="checkbox"/>	
4	<input type="checkbox"/>	Video 4	Audio 4	<input type="checkbox"/>	

Identify


Video Settings

Refresh

To start recording, press the **Start recording** button.

	Mute	Volume
	<input checked="" type="checkbox"/>	<div><div></div></div>

Video Settings Refresh



A pop-up window will appear, informing you about Noldus's Data Transparency Policy. Read it carefully before clicking the **I will inform the participants** button. You can also check the box **Don't show me again** to prevent further occurrence.

Data transparency

Responsible Use

Informed consent

You are starting the recording. Before you continue, if your study involves human subjects, please inform the participants in your study what you will be doing with their data and why. Make sure you have their consent to take part in the research.

Data Transparency

By clicking the OK button, you acknowledge your responsibility to inform the participants about:

- Which data is collected
- How the data will be used
- Where the data will be stored
- Who will have access to the data
- What the collected data will not be used for
- How the participant will benefit from the data collection

Use of MediaRecorder is subject to certain restrictions, which are documented in section 7 of the [End-User License Agreement](#).

☐ Don't show me again

I will inform the participants

The button then changes into the **Stop Recording** button and the Recording timer starts running.

2	<input checked="" type="checkbox"/>	Participant 1	Audio 3
3	<input checked="" type="checkbox"/>	Overview cam 1	Audio 2
4	<input checked="" type="checkbox"/>	Participant 2	Audio 4

Identify


Recording: 00:00:03 (hh:mm:ss)

Press the button again to stop recording. Recording stops automatically when you [specified a recording time](#) .

## how to access the recorded videos

To quickly access the recorded videos, press the **Open target folder** button in the main MediaRecorder window.




**File location:**  ... 

**File base name:**

**Output mode:**

1	2
3	



The video files are stored in the default video file location, or the location you specified. See [To select the video file location](#).

# For more information

## mediarecorder help

The MediaRecorder Help opens when you press F1 in the program. It can also be accessed in the Windows apps screen, and can be downloaded from the MyNoldus portal. From the **Help** menu select **Noldus Online** and then **MyNoldus portal**.

## training courses

Noldus offers training courses to help you get the most out of our software/hardware. In addition to on-site courses at your lab and in-house courses at our headquarters in Wageningen (The Netherlands) and Leesburg (VA, USA), we offer remote training courses.

Remote trainings are available in the same range of languages as on-site courses. Please contact your sales representative for more information about trainings.

## support

If you encounter problems, contact Noldus Technical Support, from the **Help** menu select **Noldus Online** and then **MyNoldus portal** and submit a support case or find a helpdesk in your region.