

# Connect IP cameras to the PC


## Prerequisite

You have followed the entire procedure in [Ethernet card setup](#)

## Procedure

The procedure applies to a Cisco SG 100D 8 port Gigabit PoE-switch, which is supported with MediaRecorder. Using cross-network cables (cat-5e or newer), connect the Ethernet card on the PC to one of the four right ports of the PoE-switch.

1. Connect the camera to one of the four left ports of the PoE-switch. Do not use any of the four right ports for this purpose, since they do not supply your cameras with power.
2. Power up the PoE-switch.



## Notes

- The Power over Ethernet (PoE) switch is a device that receives and sends data (in this case, video data from the camera) and passes power along on twisted-pair Ethernet cabling.
- If you use multiple IP cameras, connect them to separate ports on the PoE-switch. Make sure the capacity of the switch is high enough for the cameras.

See [Power over Ethernet](#) in [Recording devices](#)

---

Revision #6

Created 1 December 2023 08:37:47

Updated 2 January 2024 19:52:25