

How the behaviors are classified

Behavior annotation in SmartAnnotator consists of the following steps:

1. The top view videos are preliminary annotated with the Rat Behavior Recognition Module or Mouse Behavior Recognition Module of Ethovision XT version ***
See the section *Behavior Recognition Requirements* in the EthoVision XT Help for the best practices for camera setup, arena setup and detection settings.
2. SmartAnnotator makes fragments, each with a separate event, based on the automatic behavior classification in EthoVision XT.
3. Video fragments are shown to the user, for manual annotation.
4. The AI model is updated with these manual annotations.
5. Steps 3 and 4 are repeated until the classifier reaches an accuracy threshold.
6. The resulting classifier is applied to all events that are not yet annotated, giving these an automatic annotation if the model is certain enough.
7. Events that are most informative for the classifier, are shown to the user, for manual annotation.
8. Steps 3-7 are repeated, until all events are annotated.

Revision #2

Created 8 January 2025 11:20:24 by Olga Krips

Updated 9 January 2025 08:39:05 by Olga Krips