

The Viso interface

- [Log in](#)
- [Tabs](#)
- [Dashboard](#)
- [Locations](#)
- [Sessions](#)
- [Scheduler](#)
- [Codes](#)
- [Select locations for recording](#)
- [Select one location and record a session](#)
- [Select multiple locations to record separate sessions for each location](#)
- [Select multiple locations to record a composite session](#)

Log in

aim

To start working with Viso you need to log in with the login name and password that you received from your colleague who manages the Viso set-up.

procedure

1. Start Viso.
2. Enter the login details supplied by your colleague who manages the Viso set-up.

LOGIN


User name:

Password:

Login


Notes

- When you are logged in, your user name is shown in the top-right corner of the Viso window. Click it to change your password or log out.



- **TIP** Change your password the first time you log in. To do so, click your name in the top-right corner and choose **Change password**. You can only change your password if you are a local Viso user.
- Contact your colleague who manages the Viso set-up if you have lost your password.
- You can only be logged in once in Viso, independent whether this is via the regular Viso interface or via the Viso Web App. If you log in a second time, you will be logged out on the other Viso instance or webpage.

Tabs



The Viso interface has the following tabs:

- [Dashboard](#)
- [Locations](#)
- [Sessions](#)
- [Scheduler](#) (if you have the Scheduler Module)
- [Codes](#) (if you have the Annotation Module)

Viso modules and user roles

The interface depends on your Viso license and your user role. The **Scheduler** tab is present when you have the Scheduler Module, the **Codes** tab when you have the Annotation Module. The **Locations** tab is present when your user role allows session recording.

Dashboard

Viso opens with the dashboard which includes a calendar and a sessions list. Days in bold indicate when sessions that you created or received an invitation for were recorded. Click on a date in the calendar to view these sessions. A session that is being recorded or has been recorded has a link to open it.



Session color coding

Recorded sessions are marked with a color in the Dashboard and in the Scheduler. The following colors are used.

Finished

Recording

Scheduled

Force stopped

Not started/error

SESSIONS - 27 OCTOBER 2017

11:55 - 11:56

Team skills training 1 - Group study room 0.101

Mr. R. Adams

12:00 - 12:30

Team skills training 1 debriefing - Demonstration room 1

Rose Bennett

12:30 - 13:00

Team skills training 2 - Group study room 0.102

Mr. R. Adams

11:56 - 11:56

Interviewing techniques 2 --- Session stopped. The session h...

Mr. R. Adams

12:30 - 13:00

Team skills training 2 debriefing - Instruction room 0.80

Rose Bennett


Locations

The **Locations** tab is present for those whose user role allows session recording. It contains a list of the locations that are set up by your colleague who manages the Viso set-up.

DASHBOARD	LOCATIONS	SESSIONS	CODES	...	S
NAME		STATUS			
Control room		Offline			
Demonstration room		Offline			
Group study room 1.02		Locked (Rose Bennett)			
Groupb study room 0.101		Online			
Individual study room 1.76		Offline			
Instruction room 0.80		Online			
Instruction room 0.82		Recording (Rose Bennett)			
Observation room 1.20		Offline			
Waiting room		Offline			

Double-click a location to open it and start recording.

See [Sessions](#)



Location status

Locations can have the following status:

- Online – Available for recording.
- Locked – The Location is opened by another user. Administrators and Technicians can open the location, but cannot start session recording. Trainers and Trainees can only open the location if they received an invitation for the session.
- Recording – Recording takes place in that location. Administrators and Technicians can open the location, but cannot start session recording. Trainers and Trainees can only open the location if they received an invitation for the session.
- Offline and grayed out – The location has not been selected for recording by your colleague who manages the Viso set-up.
- Configuration warning – One camera, or some of the cameras in the location are not configured correctly. Consult your colleague who manages the Viso set-up to solve this. You can record with the cameras that are online.
- Configuration error and grayed out – None of the cameras in the location are configured correctly. Consult your colleague who manages the Viso set-up to solve this. You cannot open the location.

NOTE For Technicians, the Locations tab contains an Edit locations button. Technicians can add and edit locations. See the Viso Setup Manual for details.

Sessions

The **Sessions** tab contains the list with recorded sessions that you are entitled to review. For Technicians and Administrators, the **Session** tab lists all recorded sessions. For Trainers, only the sessions that they created themselves, received session ownership of, were assigned as annotator for or received an invitation for are listed. For Trainees, the list contains the sessions they received ownership of, were assigned as annotator for or received an invitation for.

Viso

MICHAEL ANDERSON

DASHBOARD

LOCATIONS

SESSIONS

SCHEDULER

CODES

Drag and drop a header here for grouping				
NAME	DATE	LOCATION	OWNER	ANNOTATOR
Room 2 - session 1	11/4/2022 10:54 AM	Room 2	Michael Anderson	Rose Bennett
Room 1 - session 2	11/4/2022 10:48 AM	Room 1	Michael Anderson	Rose Bennett
Room 1 - session 1	11/4/2022 10:47 AM	Room 1	Michael Anderson	Rose Bennett

Double-click a session to open it for review.

See [Review sessions](#)

Viso

Mok counseling session 3

BACK

EMMA SANCHEZ

INSTRUCTION ROOM 0.80

00:00:00

00:29:59

Scheduler


The Scheduler tab is available if you have the Scheduler Module. In this tab you can schedule your session recordings in advance.

The screenshot shows the Viso Scheduler interface. At the top, there's a header with the Viso logo, the user name 'ROSE BENNETT', and window controls. Below the header is a navigation bar with tabs: DASHBOARD, LOCATIONS, SESSIONS, SCHEDULER (highlighted), and CODES. The main area displays a calendar view for '06 November 2017'. The calendar has a grid with columns for different rooms: DEMONSTRATION ROOM, GROUP STUDY ROOM 0.101, GROUP STUDY ROOM 0.102, and INSTRUMENT ROOM. The rows represent time slots from 07:00 to 13:00. A session titled 'Demonstration group 1' is scheduled for 12:00 in the DEMONSTRATION ROOM, with corresponding 'Scheduled session' entries in the other rooms. The interface also includes a 'New schedule' button and a 'Today' button at the top left of the calendar area.


Day	Week	Work Week	Month	6 November 2017
	DEMONSTRATION ROOM	GROUP STUDY ROOM 0.101	GROUP STUDY ROOM 0.102	INSTRUMENT ROOM
	06 Monday	06 Monday	06 Monday	06 Monday
07:00				
08:00				
09:00				
10:00				
11:00				
12:00	Demonstration group 1	Scheduled session	Scheduled session	Scheduled session
13:00			Scheduled session	Scheduled session



Codes

The **Codes** tab is available if you have the Annotation Module. In the **Codes** tab you can create Codes lists with Subjects and Markers. This can be done by Technicians, Administrators and Trainers.

 Viso

DASHBOARDLOCATIONSSESSIONSCODES

 New list

NAME	CREATED BY
 Conversation techniques	Emma Sanchez
 Teaching techniques	Emma Sanchez

Use the codes list to annotate your sessions.

See [Annotations](#)

Select locations for recording

What do you want to do?

- [Select one location and record a session](#)
- [Select multiple locations to record separate sessions for each location](#)
- [Select multiple locations to record a composite session](#)

If you have the Scheduler Module, see also [Schedule Session Recording](#).

Select one location and record a session

aim

To select one location for a single session recording.

prerequisite


Your user role allows session recording (Technician, Administrator, or Trainer).

procedure

1. In the **Locations** tab, double-click the **name** of the location in which you want to start recording.

DASHBOARD	LOCATIONS	SESSIONS	CODES	5
NAME		STATUS		
Control room		Offline		
Demonstration room		Offline		
Group study room 1.02		Locked (Rose Bennett)		
Groupb study room 0.101		Online		
Individual study room 1.76		Offline		
Instruction room 0.80		Online		
Instruction room 0.82		Recording (Rose Bennett)		
Observation room 1.20		Offline		
Waiting room		Offline		

2. A preview of the cameras in the location appears. The toolbar displays the location name.



3. Continue with the third chapter of the Reference ([Prepare session recording](#))

NOTE The location is locked for recording by others at the moment you select it in the Locations list.

NOTE You can also record a single session with the Viso Web App, see [The Viso Web App](#)

Select multiple locations to record separate sessions for each location

aim


To select multiple locations for simultaneous session recording. Each location results in a separate session.

Prerequisites

- The setup account enabled multiple-location recording.
- Your user role allows session recording (Technician, Administrator, or Trainer).

procedure


1. In the **Locations** tab, select the locations.
2. Click the **Separate sessions** button on the toolbar.




3. A preview of the cameras in the location that you selected first appears. The toolbar displays the carousel control.

The carousel control

The location you selected first is displayed in the middle of the carousel control.



The location displayed on the right is the location you selected second.




The location displayed on the left is the location you selected last.



To open locations with the carousel control

Do one of the following:

- Press the carousel control buttons to browse through the locations.
- Press the **Page up** and **Page down** button on your keyboard.
- Click the down arrow next to the location in the middle and select another location from the list.



To set up the locations

See the third chapter of this Reference ([Prepare session recording](#)) for the procedure. Follow the procedure for each location.

NOTE If you selected only two locations, the second is displayed on the left as well as on the right in the carousel control.

NOTE The **Separate Sessions** button is enabled once you selected multiple locations. Separate session recording is supported with a maximum number of eight locations with two cameras each.

Select multiple locations to record a composite session

aim


To select multiple locations for simultaneous recording. All recordings are combined in one session.

prerequisites

- The Technician account enabled multiple-location recording.
- Your user role allows session recording (Technician, Administrator, or Trainer).

procedure

1. In the **Locations** tab, select the locations.
2. Click the **Composite session** button on the toolbar.



A preview of the cameras in the location that you selected first appears. The toolbar displays the carousel control.

See [The carousel control](#) and [To open locations with the carousel control](#) control to set up the locations

See the third chapter of this Reference ([Prepare session recording](#)) for the procedure. Follow the procedure for each location.

NOTE The Composite Session button is enabled once you selected multiple locations.

NOTE Composite session recording is supported with a maximum number of eight locations with two cameras each.