

Create and edit groups

aim


To give groups of users identical rights, for instance, permission to view your session live or to review it after it has been recorded. As a Trainer you can, for instance, create groups of students (Trainees).





prerequisites

Your user role (Technician, Administrator, Trainer) allows creating groups.

procedure

1. Click the menu button and choose **Groups**.
2. If you wish to create a new group, click **Add group** and enter a name for the group of users. If you wish to edit an existing group, click on the editing icon in its row.



GROUP	OWNER	MEMBERS	
Group 2	Hajnalka Trefas	2	 
Group 3	Hajnalka Trefas	1	 

3. Select whether the group is private or not. If you want, for instance, Trainer B to see your groups, make the groups public (leave the **Private group** checkbox unselected) and make Trainer B a group member. Trainer B can then invite your groups to view his/her sessions.
4. Select the **Department** where your group will belong.
5. From the list of users within the chosen Department, select the ones you want to add to the group, using the checkboxes. For multiple selection, use the Ctrl or the Shift key the usual way and then press the space key. It will change the status of the checkboxes of all the selected users. When you wish to remove users from a group, follow the same steps and uncheck the relevant boxes. Please note the following:
 - As an Administrator or Trainer, you will only see the users of your own department

- As an Administrator, you will not be able to add yourself to a public group; however, you can add yourself to a private group if you are already a member of the group. If you wish to be added to a group, you must be added by an Administrator. Please note, you can only be added to one group at a time. In the same way, you can only be removed from one group at a time.

6. Click **Save**.

NOTE To give other users rights, for instance, to view your session live or to review it afterwards, you have to invite them, please see [Invite participants](#).

NOTE As a Trainer, you cannot edit the public groups created by other users; however, you can use these.