

# Edit scored subjects, markers and remarks

## aim

To correct annotation errors.


## prerequisites

- Your license includes the Annotation Module.
- Your user role allows session recording (Technician, Administrator, or Trainer).
- If you have the user role Trainer, you created the session, received ownership or annotator rights. Trainees can also receive ownership or annotator rights.

## procedure


1. Open a recorded session (see [Review sessions](#)) and play it.
2. To make new annotations, enter them following the procedure in [Annotate](#)

To add a remark, press **Ctrl + R** or click this button.




To delete a scored subject, marker or remark, right-click it in the **Event Log** and click **Delete**.

To modify a scored subject or marker, right-click it and select **Modify subject**, or **Modify marker**. Then choose the correct one.



3. To save the session before you exit it, press **Ctrl + S**.
4. When done, choose **Back**. Choose **Save** to save new edits.



**NOTE** Others can view the scored markers and remarks when a recorded session is edited. To do so, open the Sessions tab and open the session. The markers and remarks appear when the person that edits the session saves it (**Ctrl + S**).

**NOTE** Scored [Markers](#) and [Remarks](#) can only be edited by the user that created the session and users that received session ownership or annotator rights.

---

Revision #3

Created 15 January 2024 11:38:49

Updated 11 March 2024 14:26:38