

# Manage scheduled sessions

## aim

To reschedule or remove a scheduled session.

## prerequisites

- Your Viso license includes the Scheduler Module.
- You created the scheduled session, or your user role is Technician, or Administrator.

## Modify a session

1. Open the **Scheduler** tab.
2. Double-click the scheduled session.
3. Edit the Schedule and click **OK**.

## Delete a session

1. Open the **Scheduler** tab.
2. Select the schedule and press the **Delete** key on your keyboard or the X icon.

**NOTE** You can only delete a scheduled session in the Scheduler. Not a session that has been recorded. To delete a recorded session, see [Manage sessions](#).

---

Revision #4

Created 15 January 2024 11:38:51

Updated 11 March 2024 14:26:38