

View annotations

aim

To view the annotations from a session recorded by someone else.

prerequisites

- Your license includes the Annotation Module.
- You either have rights to view the session (Technician, Administrator), you created the session (Trainer), you have received ownership or annotator rights or you are invited to view the session.
- The session has been recorded.

procedure

1. Open the **Sessions** tab.
2. Double-click a session to open it.
3. Play the session.

The annotations are shown in the [Event log](#) and the [Timeline](#)

Event log


To view the Event log, open the **Events** tab at the bottom of your window. To jump to an event in the Timeline and videos, click it in the Event log.

| Time | Participant | Event | Icon |
|----------|-------------|----------------|------|
| 00:01:09 | Student | Laughing | = |
| 00:01:14 | Student | Making eye... | ✱ |
| 00:01:17 | Teacher | Making eye... | ✱ |
| 00:01:25 | Student | Asking-clos... | + |
| 00:01:30 | Teacher | Informing | + |
| 00:01:48 | Student | Asking-open | ■ |
| 00:02:14 | Teacher | Answering | ■ |


The screenshot also shows a 'TIMELINE' view with a horizontal axis from 00:00:00 to 00:02:00. Markers for the events are placed along this axis, with a red box highlighting the 'Making eye...' event at 00:01:14. The interface includes tabs for 'CODES', 'EVENTS', and 'PARTICIPANTS' at the bottom.

Timeline


- To jump to an event in the videos and Event log, click it in the Timeline.
- To display an event, hover over a symbol.



- To move to the next event, click the **Skip to next event** button or press **Ctrl + right arrow** on your keyboard.



- To move to the previous event press **Ctrl + left arrow** on your keyboard.
- To zoom in and out the Timeline, hover with your mouse cursor over the Timeline and use your mouse wheel.
- To view an interval in detail, drag the borders of the line below the Timeline.



Revision #3

Created 15 January 2024 11:38:49

Updated 11 March 2024 14:26:38