Install MediaRecorder

IMPORTANT - If you have an older version of MediaRecorder - Uninstall MediaRecorder 6.0 and older versions before you install MediaRecorder 6.5.

To install MediaRecorder 6.5

1. Download the installation files from my.noldus.com, and extract the contents.
2. Double-click the file **MediaRecorder [Version number] Setup.exe**.
3. Under Installation type, select **Standard**.
4. In the **Drivers and tools** field, choose the driver of your video device.
5. Click **Next**, accept the End-User License Agreement and click Install.

Notes

- See Specifications for the supported devices.
- See Recording Devices for specific information on your device like the maximum number of cameras you can use simultaneously and the maximum supported frame rate and resolution.

Revision #11
Created 2023-11-29 14:23:05 UTC
Updated 2023-12-20 16:41:29 UTC