Camera setup

Aim

To configure the cameras for use with MediaRecorder.

Prerequisites

- You have Axis IP cameras that are supported with MediaRecorder.

  See Supported IP cameras in Recording devices
- You followed the procedure in Ethernet card setup
- You connected all cameras.

Procedure

You must carry out the following steps for each camera:

1. Set camera IP address
2. Set camera password and select power line frequency
3. Create camera settings
4. Create audio settings

Revision #5
Created 2023-12-01 08:38:50 UTC
Updated 2023-12-20 12:47:08 UTC