Create The Observer XT settings

Aim


To enable communication between The Observer XT and N-Linx Server.

Prerequisites

- You installed N-Linx Server on a computer in the network.
- If MediaRecorder and The Observer XT, or N-Linx Server run on different computers, you created exceptions for the port that N-Linx uses in Windows Firewall.
  See Create exceptions for N-Linx port in Windows Firewall
- N-Linx Agent is installed on the computer with MediaRecorder and connected to N-Linx Server.
  See Install N-Linx Agent
- MediaRecorder is open.


Procedure

1. In The Observer XT, choose **File** > **Preferences** > N-**Linx Settings**.
2. Select **Use N-Linx server to connect with other applications.**




3. In the N**-Linx Server address** field, enter the computer name of the computer with N-Linx server. If The Observer XT is installed on the same computer as N-Linx Server, leave the default entry Localhost.
4. In the **N-Linx Server port** field, leave the default port 5672. Ask your system administrator for assistance if this port is used by another program.
5. Click **Test connection**. If connection with N-Linx Server is found, the **Status** will change to **Connected**.
6. Click **OK**.


7. Choose **Setup** > **Project**
8. Under Observation source not, click the **Devices** button.



9. Select **MediaRecorder [Version number][Computer name]** and click **Edit Settings**. MediaRecorder now opens if it is not already open.



0. Click **Edit Settings** again. A window with the available configurations (see Save the camera settings in Record video) now opens. Select the configuration file you want to use and click **OK**. If you do not select a configuration, **Current Settings.mrs** will be used. Optionally enter a unique name in the **Identify device** as field. This way you can control multiple instances of MediaRecorder with one instance of The Observer XT.



1. Repeat the steps 9 and 10 for each MediaRecorder instance.

Note

When you close MediaRecorder, the MediaRecorder line in the **Devices** window (step 9) will be grayed out and the status will change to **Disabled**. MediaRecorder will start when you create an observation.

Revision #9
Created 2023-11-29 22:20:56 UTC
Updated 2023-12-20 12:47:09 UTC