The steps to install MediaRecorder

IMPORTANT

If you have an older version of MediaRecorder - Uninstall MediaRecorder 6.0 and older versions before you install MediaRecorder 6.5.

If you have a hardware key - Do not insert your key before the installation is complete. If you do so, other drivers are automatically installed, which may lead to incorrect functioning of the license key. Drivers of the license key should be installed using the installation USB stick.

See Install MediaRecorder and device drivers

Procedure

1. Turn off automatic updates for device drivers
2. Select correct power options
3. Install MediaRecorder and device drivers

Revision #8
Created 2023-11-29 16:18:13 UTC
Updated 2023-12-20 21:21:58 UTC