Specify a default session name

aim


To set a default name for your sessions. Your sessions will be saved with the default name, followed by a consecutive number.

prerequisite


- Your user role allows session recording (Technician, Administrator, or Trainer).

procedure

1. Click the **Menu button** at the top of the Viso window and choose **Settings > Default Session name**.



2. Select **Save session with a default name** and enter a name in the **Default name** field.



3. Optionally, select **Edit default name before the session** is saved. This way you can give the session another name than the default one with a consecutive number. This is, for example, useful if you want to add the name of the subject to the session name.

Revision #2
Created 2024-01-15 10:38:45 UTC
Updated 2024-03-08 10:02:28 UTC