

# System requirements

## Operating system

MediaRecorder 6 supports Windows 11 64 bit Professional edition.

## Language packs

MediaRecorder 6.5 was thoroughly tested with the Windows 11 U.S. English language packs. It is possible that certain local language versions of Windows may affect the performance of the program.

## Touch features

MediaRecorder is not designed for use with Windows touch features.

## Computer

### Desktop computer

MediaRecorder 6.5 has been tested with a Dell Precision™ T3640 and 5820 PC. If you order a complete solution from Noldus Information Technology, you will obtain one of these computers, or a successor, with MediaRecorder installed and ready to use.

### Laptop

MediaRecorder 6.5 has been tested on a Dell Precision™ 3551 mobile workstation with the following devices only:

- 2 Axis P1375 cameras (Portable Observation Lab). For more information on the Portable Observation Lab, consult the Reference Manual on the desktop of the notebook.
- 1 Axis P1375 camera and 1 Epiphan Pearl Nano (Portable Usability Lab). For more information on the Portable Usability Lab, consult the Reference Manual on the desktop of the notebook.
- 1 Basler acA1920-155um USB 3.0 camera.

## Technical specifications test computers

Dell Precision™ T3640 PC

- Processor: Intel i7-10700, 4.8 GHz.
- Internal memory: 8 GB.

- Hard disk: 1 TB.
- Graphics card: 5 GB NVIDIA Quadro P2200.

#### Dell Precision 5820

- Processor: Intel i9-10900X CPU @3.70 GHz
- Installed RAM: 8.00 GB
- Graphics card: Nvidia Quadro Pro P2200

#### Dell Precision™ 3551 mobile workstation

- Processor: I7-10750H (6-core), 2.2 GHz.
- Internal memory: 8 GB.
- Graphics card: 4 GB Nvidia Quadro P620 w.
- Hard disk: 1TB.

---

Revision #5

Created 29 November 2023 16:50:41

Updated 21 December 2023 10:16:46 by Niek Wilmink