

# What's new compared to version 6.2

## Selection of audio source for all the videos

While previously only Audio device 1 could be used as a common audio source for all the videos, now in Picture-in-Picture mode, you may select which audio source to use for this purpose.

## Data Transparency Policy

When starting recording, a pop-up window will appear, informing you about Noldus's Data Transparency Policy.

If you have been using a previous version of MediaRecorder, you should read these topics to get an idea of the improvements in this version.

- [What's new compared to 6.0](#)
- [What's new compared to version 5](#)

---

Revision #8

Created 29 November 2023 15:28:18

Updated 20 December 2023 21:01:26