


Codes

The **Codes** tab is available if you have the Annotation Module. In the **Codes** tab you can create Codes lists with Subjects and Markers. This can be done by Technicians, Administrators and Trainers.



Viso

DASHBOARD LOCATIONS SESSIONS **CODES**

New list

NAME	CREATED BY
Conversation techniques	Emma Sanchez
Teaching techniques	Emma Sanchez

Use the codes list to annotate your sessions.

See [Annotations](#)

Revision #3

Created 15 January 2024 11:38:46

Updated 11 March 2024 14:26:38