

# Dashboard

Viso opens with the dashboard which includes a calendar and a sessions list. Days in bold indicate when sessions that you created or received an invitation for were recorded. Click on a date in the calendar to view these sessions. A session that is being recorded or has been recorded has a link to open it.



## Session color coding

Recorded sessions are marked with a color in the Dashboard and in the Scheduler. The following colors are used.

	SESSIONS - 27 OCTOBER 2017
Finished	<div><div></div><div>11:55 - 11:56</div><div>Team skills training 1 - Group study room 0.101</div><div>Mr. R. Adams</div></div>
Recording	<div><div></div><div>12:00 - 12:30</div><div>Team skills training 1 debriefing - Demonstration room 1</div><div>Rose Bennett</div></div>
Scheduled	<div><div></div><div>12:30 - 13:00</div><div>Team skills training 2 - Group study room 0.102</div><div>Mr. R. Adams</div></div>
Force stopped	<div><div></div><div>11:56 - 11:56</div><div>Interviewing techniques 2 --- Session stopped. The session h...</div><div>Mr. R. Adams</div></div>
Not started/error	<div><div></div><div>12:30 - 13:00</div><div>Team skills training 2 debriefing - Instruction room 0.80</div><div>Rose Bennett</div></div>