

# Sessions

The **Sessions** tab contains the list with recorded sessions that you are entitled to review. For Technicians and Administrators, the **Session** tab lists all recorded sessions. For Trainers, only the sessions that they created themselves, received session ownership of, were assigned as annotator for or received an invitation for are listed. For Trainees, the list contains the sessions they received ownership of, were assigned as annotator for or received an invitation for.

Viso

MICHAEL ANDERSON

DASHBOARD

LOCATIONS

SESSIONS

SCHEDULER

CODES

Drag and drop a header here for grouping				
NAME	DATE	LOCATION	OWNER	ANNOTATOR
Room 2 - session 1	11/4/2022 10:54 AM	Room 2	Michael Anderson	Rose Bennett
Room 1 - session 2	11/4/2022 10:48 AM	Room 1	Michael Anderson	Rose Bennett
Room 1 - session 1	11/4/2022 10:47 AM	Room 1	Michael Anderson	Rose Bennett

Double-click a session to open it for review.

See [Review sessions](#)


Viso

Mok counseling session 3

BACK

EMMA SANCHEZ

INSTRUCTION ROOM 0.80



00:00:0000:29:59

