


# Tabs



The Viso interface has the following tabs:

- [Dashboard](#)
- [Locations](#)
- [Sessions](#)
- [Scheduler](#) (if you have the Scheduler Module)
- [Codes](#) (if you have the Annotation Module)

## Viso modules and user roles

The interface depends on your Viso license and your user role. The **Scheduler** tab is present when you have the Scheduler Module, the **Codes** tab when you have the Annotation Module. The **Locations** tab is present when your user role allows session recording.

Revision #3

Created 15 January 2024 10:38:46

Updated 11 March 2024 13:26:38