

# What's new in Viso 12?

## ***Software based licensing fully integrated in Viso***

You can activate or deactivate your software license from the Viso application directly.

The license key needs to be able to check its validity online at least every 90 days. The 90 days check counter is reset after a successful license check. If the license check has not been successful after 90 days, your license becomes suspended. Please contact support when this happens.

## ***Highlight videos***

In the Annotation module, you can now create highlight videos which summarise the most interesting point events of your sessions to be used e.g. in presentations or lectures. To these highlight videos, you may add your markers and remarks as subtitles.

## ***The right to manage recording options moved to Technicians***

Those with a Technician account can now enable multiple location recording, including composite and separate sessions.

## ***Allocate groups to departments***

When creating groups, now the department where the group belongs can also be provided.

## ***Data Transparency Policy***

When a new user logs in, a pop-up window will appear, informing the user about Noldus's Data Transparency Policy.

---

Revision #4

Created 4 March 2024 08:41:19

Updated 21 January 2025 17:17:36