What’s new compared to version 6.2

Selection of audio source for all the videos

While previously only Audio device 1 could be used as a common audio source for all the videos, now in Picture-in-Picture mode, you may select which audio source to use for this purpose.

Data Transparency Policy

When starting recording, a pop-up window will appear, informing you about Noldus's Data Transparency Policy.

If you have been using a previous version of MediaRecorder, you should read these topics to get an idea of the improvements in this version.

- What's new compared to 6.0
- What’s new compared to version 5

Revision #8
Created 2023-11-29 15:28:18 UTC
Updated 2023-12-20 21:01:26 UTC