
How Viso annotations appear in The Observer XT

In the Coding Scheme

In The Observer XT, choose **Setup > Open Coding Scheme**, to open the Coding Scheme. Viso subjects appear in the **Subjects** panel. Viso markers are added as Point events in a Start-Stop group. Remarks are not added to the Coding Scheme.



In the observation

Choose **Observe > Observation > Open > [session name and time]** to view the annotations in the Event Log. Viso Subjects are also Subjects in The Observer XT. Markers are point events and are indicated with blue dots. Remarks are imported as comments in the Event Log.

Time	Subject	Behavior	Comment
0.00	**Start**		
12.06	Student	Talking	
17.46	Teacher	Asking-open	
24.34	Student	Listening	
30.62	Student	Answering	
39.22	Teacher	Instructing	
49.82	Student	Laughing	
55.75	Teacher	Laughing	
59.31	Student	Asking-closed	
67.75	Teacher	Summarizing	
75.03	Student	Making eye cont	
80.47	Teacher	Making eye cont	
86.27			Disturbance in back of classroom
99.03	Teacher	Checking	
106.08	Student	Summarizing	
108.24	Teacher	Asking-open	
116.44	Student	Answering	
125.74	**Stop**		

To jump to the video frame at which a certain marker or remark was scored, click its row in the Event Log.

Revision #1
Created 2024-01-15 10:38:52 UTC
Updated 2024-01-15 10:43:02 UTC