
Tabs



The Viso interface has the following tabs:

- Dashboard
- Locations
- Sessions
- Scheduler (if you have the Scheduler Module)
- Codes (if you have the Annotation Module)

Viso modules and user roles

The interface depends on your Viso license and your user role. The **Scheduler** tab is present when you have the Scheduler Module, the **Codes** tab when you have the Annotation Module. The **Locations** tab is present when your user role allows session recording.

Revision #3
Created 2024-01-15 10:38:46 UTC
Updated 2024-03-11 13:26:38 UTC